


PowerTrack



All in one pressure sensors treadmill, removable sensors for walkway use.

PowerTrack is the result of 37 years of developments in the field of baro-podometry. Unique on the market, PowerTrack mix our professional pressure sensors to a strong gym treadmill. Powertrack allows various type of assessments: from fall risk patient to top athlete.







4 types of acquisitions



- Static
- Ground Dynamic
- Posturologic
- Treadmill Dynamic

Software release 10 , Unmatched features :

- Thermographic view , isopressure, 3D, histograms, diagrams,
- Exclusive features : FPI (Foot Posture Index), CPEI (Center of Pressure Excursion Index), Diabetic foot pressure Index , Ulcerations level of risk by David Armstrong.
- Sway index , COP Butterfly, symetry, balance, pace, angles etc...
- Customizable clinical assessment.
- Customizable reports, PDF export.
- Database importation of previous Medicapteurs software
- Vidéo synchronisation with to 2 cameras.



Size (length/height/width)	180 x153 x 88 cm
weight	125 kg
Speed	0,8 to 20 km/h
Pressure min/max per sensor	0,6 N to 100 N
PC connection	USB (4.5 meter wire)
Acquisition frequency	Up 200 images/second